

## Clinitemp™

---

### Product Data Sheet

---

**Product Description:** L314 Latex Free Forehead Moving Line Temperature Monitor

---

**Standard Range:** 84-106°F / 29-41°C

---

#### Construction:

Adhesive Tape	Adhesive: Acrylic copolymer pressure sensitive adhesive, thickness 37μ each side of carrier Carrier: Non-woven tissue, thickness 33μ Release Liner: Fully bleached, thickness 178μ Operating temp: -1 to 82°C (continuous use)
Covering Film	Polyester
Color Change Material	Nontoxic mixture of liquid crystal esters
Accuracy	± 1.00°F / 0.6°C
Shelf Life	12 months from invoice date when stored at room temperature & humidity (i.e. 21°C (70°F) & 50% relative humidity).

---

**Application Data:** Peel temperature-indicating label from backing paper. Apply to forehead. Reversible analog line moves along the scale to show detailed temperature readings in °F and °C.

**Function:** Continuous display of temperature, as indicated by a bright line, of anesthetized patients. Core adjusted with 3.5°F / 2°C offset.

**Use:** Ideal for use in Pre-Op, Surgery, PACU